

EZee Life™

16" Bathtub Safety Rail

Model: CH2034



EZee Life™ Products

34 Futurity Gate, Unit 15
Concord, Ontario L4K 1S6



Limited Warranty:

EZee Life™ products carry a one year warranty against manufacturing defects, faulty materials and workmanship only.

Any alterations, misuse, abuse, or accidental damage voids this warranty. EZee Life™ reserves the right to replace or repair any part that has become defective.

The warranty does not extend to consequential costs resulting from the fault or defect of the product: freight and travel costs, loss of earnings, or other expenses that one may incur.

Warranty is valid for the original purchaser only and the original receipt must be presented before any warranty options are considered. Please also be prepared to provide the date of purchase and serial number.

Specifications:


Height (when mounted):	16"
Secondary grip height:	9"
Bar diameter:	1.75"
Installation method:	Clamp
Clamp opening:	2 3/4" to 6 1/2"
Product weight:	5.1 lbs 2.32 kgs
Weight capacity:	300lbs (may vary depending on clamp tightness)

Features:

- Easy to clean plastic
- Ergonomic design
- Safe grip
- Easy to install no tools needed
- Won't rust or scratch (clamp excluded)
- Attractive design
- Dual grip design

Installation:

1. Clean and dry the area of the bathtub wall where the bar is to be installed to ensure a secure grip.
2. Fold out clamp handle and unscrew to widen clamp slightly wider than tub wall.
3. Lower the tub bar all the way down onto tub wall in the desired position with the clamp handle positioned outside of the tub.
4. Tighten the clamp using the clamp handle until sufficiently tight enough to prevent the bar from slipping during use.
5. Test the stability of the bar by pulling forward, backwards and side to side. If the bar slips any amount ensure the tub wall is not slippery due to soap film or other substance and re-tighten the clamp.



Tub safety rails normally require retightening from time to time. Check periodically before use and re-tighten as needed.

Over-tightening the tub rail can lead to damage to the safety rail or the tub wall.

CAUTION: Tub safety rails are not suitable for all bath tub styles. Tubs with excessively sloped walls or top lips may not allow the bar to be securely installed. Fibreglass tubs may not have solid enough walls to withstand the pressure of the clamp.

Directions for Use:

The upper grip area of the tub rail is intended to assist the user's balance when stepping into and out of the tub. It is not intended to support the user's full weight at any time.

The lower grip area is more secure and can be used to support the user's weight when rising out of or lowering into the bath tub.



Safety Information:

DO NOT install this equipment without first reading and understanding this user manual. If you are unable to understand the warnings, cautions and instructions contact your dealer.

Hold the tub bar securely when entering or exiting the bath tub.

DO NOT use this product if you exceed its weight capacity.

ALWAYS ensure that the tub rail is secure prior to use.

DO NOT use this product if there is any apparent damage.

DO NOT exert full weight on upper grip area.

Cleaning Information:

EZee Life products can be cleaned with most commercial cleaning products. Do not use abrasives.

We recommend cleaning the surfaces of the product and attachments with a solution made from a mild detergent (dish soap, etc.) to remove organic materials.

Once cleaned, it can be disinfected by wiping the product down with a 5% (1 cup bleach to 1 gallon water) solution of bleach and water. Allow the bleach solution to remain on the surfaces for 3 to 5 minutes, then rinse thoroughly with clear water.

Note: Bleach reacts with ammonia and organic matter (feces) causing gases that can irritate the eyes and lungs. Use in a ventilated area and avoid contact with other chemicals and organic matter.

Imported / Distributed by:

Factory Direct Medical
34 Futurity Gate, Unit # 15
Concord Ontario L4K 1S6
Local: 416-739-8393
Toll Free: 855-235-2400