

EZee Life™

Folding Rollator

Models: CH3016 (Black)



EZee Life™ Products

15-34 Futurity Gate, Concord, Ontario L4K 1S6



Limited Warranty:

Ezee Life™ products carry a one year warranty against manufacturing defects, faulty materials and workmanship only. Any alterations, misuse, abuse, or accidental damage voids this warranty. Ezee Life™ reserves the right to replace or repair any part that has become defective.

The warranty does not extend to consequential costs resulting from the fault or defect of the product: freight and travel costs, loss of earnings, or other expenses that one may incur.

Warranty is valid for the original purchaser only and the original receipt must be presented before any warranty options are considered. Please also be prepared to provide the date of purchase and serial number.

General Safety Warning:

Please ensure these instructions are fully read and implemented. Failure to do so may result in injury to the user. Please keep these instructions in a safe location for future use.

Intended Use:

The lightweight 4 wheel rollator is designed to provide support and assistance to those unable to walk with some level of support. This rollator is designed for both indoor and outdoor use and can be folded for easy storage and transportation.

Safety Precautions:

- Do not Exceed a maximum user weight of 136 kg (300 lb)
- Do not use on stairs or escalators
- Do not push the rollator with someone sitting on the seat
- Do not use the seat to transport people or objects
- Do not use the seat without the backrest in place

Safety Precautions (con't):

- Do not sit on the seat while the rollator is on an incline
- Only use the rollator as a walking support
- Only use specified part and accessories for this rollator
- Loads should only be transported in the basket
- Do not hang heavy items from the push handles

Safety Precautions - Pre Use Checks

- Ensure rollator is locked open before use
- Check the brake function before use
- Ensure moving parts are secure and working properly

Safety Precautions – During Use

- Do not self-propel the rollator while sitting on the seat.
- Serious injury to the user and damage to the rollator may result from improper use
- Do not use rollator to walk backwards
- Use caution when negotiating curbs, slopes, slippery and uneven surfaces
- Do not do any adjustments to the rollator while in use
- Use caution when the basket is loaded (5 kg / 11 lb max.)
- Brakes must be locked before using seat
- Brakes must remain locked while the user is seated





Contents List:

Item	Description	Qty.
1	Rollator Frame	1
2	Storage Basket (fitted to frame)	1
3	Pentagon Hand Bolt	2
4	Backrest	1
5	Front Wheels	2
6	Rear Wheels with Brake & Cable	2
7	Drink Holder & Extended Machine Screw	1
8	User Instruction Manual	1

Unpacking:

- Carefully check all parts prior to assembly
- Do not assemble if parts are damaged or missing
- Ensure handles and brake cables are not twisted to prevent brake problems

Assembly & Height Adjustment:


NOTE: Before assembly consult a suitably qualified healthcare professional regarding the correct height adjustment for maximum support and correct brake activation.

1. To assemble; insert the front & rear wheels using the button clips on the rear wheels and the “C” clips on the front wheels. Loosely insert the front handlebars (do not fix the handlebars at this stage).
2. Fit the backrest into the frame and secure with snap clips.
3. To set the seat height; with the rollator frame open, push down on the set until all four wheels are extended and the cross member is in line s shown in diagram (1). The seat height can be adjusted by sliding the wheel assemblies in the frame tubes and positioning using the button & “C” clips. Make sure that all the clips “click” into position and all four wheels make contact with the floor

Assembly & Height Adjustment (con't):

4. Set the correct height of the handle bars as follows:
 - Stand behind the rollator with your shoulders relaxed and your arms by your sides, looking straight ahead
 - In this position, the handgrips or the rollator should be adjusted to the height level of your wrists. (see diagram 1)
 - This will position your arms with approximately 20° to 30° bend when using the rollator
 - Check that the holes in the handle bars and the frame are aligned.
 - Place supplied bolts through the aligned holes in the handlebar and frame then hand tighten firmly but do not over tighten.
5. Check the brakes thoroughly before use. NOTE: The brakes are factory set and if assembled correctly do not need adjusting prior to use. To check; Apply and release the parking brake on each side ten times (push the brake lever down to lock, lift to release). The brake should grip the wheel and prevent any movement when applied. When released, the gap between the brake moulding and the tire should be 3 to 4 mm. If either of the above conditions is not correct, adjust the brakes as described in the maintenance section and recheck.
6. Check all other nuts, bolts and screws are tightened.
7. Check all the wheels are in good working order.
8. Check the basket is securely fixed and all the fixing buttons are tightly fastened
9. Install drink holder:
 - 1) Remove the bolt holding the rear leg of the walker to the frame/backrest bracket
 - 2) Replace the bolt with the supplied longer bolt with the head of the bolt to the inside of the frame
 - 3) Fit the drink holder onto the threaded end of the bolt and thread the nut onto the threads and tighten





Opening the Rollator

- Hold one handlebar and push down on the seat until all four wheels are fully extended and the seat is flat.
- Before using the rollator make sure that all four wheels are making contact with the floor and the rollator is stable in the open position as in diagram 1.

Folding the Rollator

- Lift the seat and pull up the seat support bar until the rollator closes as shown in diagram 2.
- For more compact storage, pull up the backrest



Walking with the Rollator

The unit is to be used as a mobility aid, creating a natural gait as the user walks securely with their primary weight distribution between the handlebars and with the rollator in front of the user.

Sitting on the Rollator

When using the rollator as a seat, ensure the safety precautions detailed on page one are adhered to. In particular:

- Ensure the parking brakes are locked
- DO NOT sit on the seat unless the backrest is in position
- DO NOT use the seat when parked on sloping or uneven ground
- DO NOT exceed maximum user weight - 136 kg (300 lb)
- DO NOT use the rollator as a wheelchair or transport device

Checking and Operating the Brakes

- Before using the rollator always check the brakes are operating properly.
- Check the normal brake function by pulling up the brake handles as shown in diagram 3a. The rollator should not roll.
- To apply the parking brakes, push down the handles until they lock down. A distinct click indicates that the brake is in the correct position. The rollator should hold firm with no movement (see diagram 3b). It is important both brakes are used for parking. Lift up on the handles to release brakes.




Diagram 3a



Diagram 3b



Brake Operation

- To slow down, pull up both brake handles (diagram 3a).
- To apply the parking brakes, push down both brake handles until they lock. A distinct "click" will be heard when the brakes are in the correct position (diagram 3b). It is important that both brakes are put in the parking position when they are used.
- Lift both handles to release the parking brakes.
- If you apply the parking brakes to stop while moving, and the rear wheels on the rollator still move, stop use immediately. Do not use again until the the brakes have been checked and correctly adjusted. Refer to the maintenance section.

Cleaning Instructions

- Clean regularly by wiping down with a soft non-abrasive cloth.
- Ensure there is no buildup of dirt or debris around the axle, brakes and wheels.

Maintenance

- Regularly check for wear and tear and ensure all fasteners are securely fastened.
- Check the brakes are working satisfactorily by both walking and squeezing the brakes and by braking to a stop.
- Check the parking brakes by pushing the handles into the locked position and then trying to push the rollator.
- If the brakes do not completely lock the wheels, the brakes need to be adjusted.

Brake Adjustment (See Diagram 4)

- There are two brake adjusters on each side of the rollator, one in front of the brake handle and the second is positioned on the rear legs behind the wheels.
- To adjust the handle brake, loosen the lock nut (you may need a tool to carry out this operation), turn the adjuster anti-clockwise to tighten the brake then re-tighten the lock nut (diagram 4).
- To adjust the rear legs brake, loosen the lock nut (you may need a tool to carry out this operation), turn the adjuster counter-clockwise to tighten the brake then retighten the lock nut. For correct adjustment the brake mechanism, should not touch the wheel but the gap should be no greater than 3 to 4 mm.
- Caution: After making any adjustments to the brakes, you must test them before the rollator is used.
- If adjustment does not tighten the brakes, do not use the rollator and contact your dealer.




Diagram: 4 – Handle brake adjustment

Imported / Distributed by: Factory Direct Medical
34 Futurity Gate, Unit # 15
Concord Ontario L4K 1S6
Local: 416-739-8393
Toll Free: 855-235-2400