

# ↑ EZee Life™ ELECTRIC CONVENIENCE VEHICLES

## PRODUCT DESCRIPTION, FEATURES & SPECIFICATIONS



# Electric Convenience Vehicles

ECVs are electrically operated mobility devices that are loaned, or rented, to patrons for short term use at Amusement Parks, Zoos, Exhibitions, National & Provincial Parks, or anywhere a lot of walking is involved and wants to be more accessible to people with mobility issues.

## Benefits to the facility of offering ECVs:

- 1) Provide facility accessibility to attract more visitors.
- 2) Create a new income stream for facilities over and above what can be generated with manual wheelchairs that often require caregiver assistance.

## Benefits to the patrons:

- 1) ECVs provide freedom and independence while visiting your facility.
- 2) Patrons can use mobility devices that are chosen to work in the terrain of the facility.
- 3) Patrons don't have to bring their own mobility devices to enjoy your facility.
- 4) Patrons can have mobility devices that may be more appropriate for the facility terrain than their own.

# C O N T E N T S

---

## Power Wheelchair ECVs

<b>Generation 1:</b> .....	4
CH4050 & CH4055	
CH4060 & CH4065	
<b>Generation 2:</b> .....	8
CH4051 & CH4056	
CH4053 & CH4054	
CH4058 & CH4059	
<b>Generation 3:</b> .....	11
CH4080	
CH4085	
<b>Generation 4:</b> .....	13
CH4075	
<b>Generation 5:</b> .....	14
CH4070	

## Mobility Scooter ECVs

<b>Travel Mobility Scooter:</b> .....	15
<b>Classic Mobility Scooter:</b> .....	16
<b>Elite Mobility Scooter:</b> .....	17
<b>EZ5S Mobility Scooter:</b> .....	18
<b>EZ6S Mobility Scooter:</b> .....	19
<b>EZEE ROVER 4 Scooter:</b> .....	20
<b>EZee Auto Fold Scooter:</b> .....	21

# G1 POWER WHEELCHAIRS

## CH4050 EZEE FOLD STANDARD

The EZee Fold Standard power wheelchair is a lightweight portable power wheelchair that folds easily for easy transport and storage but offers a range of up to 15 kms per battery. Optional second battery can be purchased to double the range.

With a seat width of 16” and a seat depth of 17” and a weight capacity of up to 264 lbs. this chair will fit most average sized adults.



SPECIFICATIONS	
Seat Width:	400 mm (16")
Seat Height:	460 mm (18.1")
Seat Depth	430 mm (17")
Backrest Height:	450 mm (17 3/4")
Net Weight:	23.1 kg (51.4 lbs.) 1 battery
Weight Capacity:	120kgs (264 lbs)
Caster:	8" PU Tire
Rear Wheel:	10" PU Tires
Battery:	Lithium battery, 24 V 9 Ah
Max speed:	6.5 kph (4 mph)
Max. Climbing Grade:	10 degrees
Charging period:	5 hours per battery
Endurance:	15 km (9.4 miles) km/battery
Overall size:	870 x 560 x 910 mm (34.2" x 22" x 35.8")
Folded size:	560 x 300 x 660 mm (22" x 11.8" x 25.9")

## CH4055 (8") EZEE FOLD

The EZee Life™ Ezee Fold CH4055 is similar to the EZee fold standard but has standard footrests. It also features a high weight capacity at 120kg and adjustable height back. This Ezee Fold model has a two battery capacity so you can double you range with the purchase of an optional second battery. With a 18.1" seat height and a 16" seat width the CH4055 will fit most tall users of average weight. Like all Ezee Fold models the CH4055 will fold in seconds and be ready for transport.



SPECIFICATIONS	
Seat Width:	400 mm (16")
Seat Height:	460 mm (18.1")
Seat Depth	430 mm (17")
Backrest Height:	450mm (17 3/4")
Net Weight:	23.1 kg (51.4 lbs.)
Weight Capacity:	120 kg (264 lbs.)
Caster:	8" PU Tire
Rear Wheel:	10" PU Tire
Battery:	Lithium battery, 24 V 9 Ah
Max speed:	6.5 kph (4 mph)
Max. Climbing Grade:	10 degrees
Charging period:	5 hours per battery
Endurance:	Up to 15 km (9.4 miles) /battery
Overall size:	870 x 560 x 910 mm (34.2" x 22" x 35.8")
Folded size:	560 x 300 x 660 mm (22" x 11.8" x 25.9" )

# G1 POWER WHEELCHAIRS

## CH4060 EZEE FOLD HD

The EZee Fold HD power wheelchair is a lightweight portable power wheelchair that folds easily for easy transport and storage but offers a range of up to 15 kms per battery. Optional second and third batteries will increase range accordingly.

With a seat width of 18“ and a seat depth of 17 and a weight capacity of up to 350 lbs. this chair will fit most larger sized adults.



SPECIFICATIONS	
Seat Width:	457 mm (18")
Seat Height:	422 mm (16.6")
Seat Depth	400mm (15 3/4")
Height adjustable backrest:	450 to 650 mm range (17.7" to 25.6")
Net Weight:	25.3 kg (56 lbs.) with 1 battery
Weight Capacity:	160 kg (350 lbs.)
Caster:	8" PU Tire
Rear Wheel:	12" Rubber/Pneumatic tire
Battery:	Lithium battery, 24 V 10 AH
Max speed:	7.2 kph (4.5 mph)
Max. Climbing Grade:	12 degrees
Charging period:	5 hours per battery
Endurance:	15 km (9.4 miles) per battery
Overall size:	1000 x 640 x 890 mm ( 39.4" x 25.2" x 35")
Folded size:	640 x 335 x 740 mm (25.2" x 13.2" x 29.1 ")

## CH4065 (12") EZEE FOLD

The EZee Life™ Ezee Fold CH4065 is similar to the Ezee fold standard but has the choice of elevating or standard footrests and has 12" drive wheels for rougher terrain. It also features a high weight capacity at 165kg and adjustable height back.

This Ezee Fold model has a two battery capacity so you can double your range with the purchase of an optional second battery. With a 20" seat height and a 17" seat width the CH4065 will fit most tall users of average weight. It also features a weight capacity of 160kg and adjustable height back.



SPECIFICATIONS	
Seat Width:	430 mm (17")
Seat Height:	507 mm (20")
Seat Depth:	430 mm (17")
Backrest Height:	441 mm (17.4")
Net Weight:	26.6 kg (58.6 lbs.)
Weight Capacity:	160 kg (350 lbs.)
Caster:	8" PU Tire
Rear Wheel:	12" Rubber Tire
Battery:	Lithium battery, 24 Volt 10 AH (1 or 2)
Max speed:	7.2 kph (4.5 mph)
Max. Climbing Grade:	12 degrees
Charging period:	5 hours per battery
Endurance:	Up to 15 km (9.4 miles) /battery
Overall size:	954 x 629 x 948 mm (37.6" x 24.8" x 37.3")
Folded size:	629 x 346 x 752 mm (24.8" x 13.6" x 29.6")

# G2 POWER WHEELCHAIRS

## CH4051 (BLACK FRAME) & CH4056 (RED FRAME) EZEE FOLD

The EZee Fold 2G electric wheelchair with 8” wheels is our 2nd Generation lightweight portable power wheelchair that folds easily for easy transport and storage but offers a range of up to 16 kms per battery. The option of being able to add a second battery will double the range.

With a seat size of 16“ x 16” and a weight capacity of up to 120 kgs (264 lbs.) this chair will fit most average sized adults.



SPECIFICATIONS	
Seat Width:	406mm (16")
Seat Height:	495mm (19.5")
Seat Depth	406mm (16")
Backrest Height:	457mm (18")
Net Weight:	52.35 lbs (23.8 kg)
Weight Capacity:	120kg (264 lb)
Caster:	178mm (7")
Rear Wheel:	203mm (8")
Battery:	Li-ion Batteries
Max speed:	6 km/h (3.75 mph)
Max. Climbing Grade:	12 degrees
Charging period:	6 - 8 Hours per Battery
Endurance:	Up to 12 km per battery
Overall size:	984 x 507 x 914 mm (38.75" x 23.5" x 36")
Folded size:	231 x 597 x 806 mm (15" x 23.5" x 31.75")

## G2 POWER WHEELCHAIRS

### CH4053 (BLACK FRAME) & CH4054 (SILVER FRAME) EZEE FOLD

This EZee Fold 2G electric wheelchair with 8" wheels is a lightweight portable power wheelchair with a folding frame for easy transport and travel. It provides a range of up to 20 km. The option of being able to add a second battery will double the range.



SPECIFICATIONS	
Seat Width:	400mm (15.75")
Seat Height:	470mm (18.5")
Seat Depth:	400mm (15.75")
Backrest Height:	495mm (19.5")
Net Weight:	21 kg (46 lbs) (single battery)
Weight Capacity:	120kg (265lbs)
Caster:	178mm (7")
Rear Wheel:	203mm (8")
Battery:	Li-ion Batteries
Max speed:	6 km/h (3.75 mph)
Max. Climbing Grade:	12 degrees
Charging period:	6 - 8 Hours per Battery
Endurance:	Up to 12 km per battery
Overall size:	889 x 597 x 914 mm (35" x 24" x 37")
Folded size:	787 x 597 x 330 mm (31" x 24" x 13")

# G2 POWER WHEELCHAIRS

## CH4058 (BLACK FRAME) & CH4059 (ORANGE FRAME) EZEE FOLD

The EZee Fold 2G electric wheelchair with 12” wheels is one of our newest 2nd Generation lightweight portable power wheelchairs with an easy folding frame for easy transport and travel that offers a range of up to 22 km.



SPECIFICATIONS	
Seat Width:	470 mm (18.5")
Seat Height (back):	470 mm (18.5")
Seat Depth	400 mm (15 3/4")
Backrest Height:	483mm (19")
Net Weight:	23 kg (50.7 lb)
Weight Capacity:	150KG (330 lbs)
Caster:	203 mm (8")
Rear Wheel:	318 mm (12.5")
Battery:	24V/7 AH x 2 PCS
Max speed:	6 km/h (3.75 mph)
Max. Climbing Grade:	12 degrees
Charging period:	6 - 8 Hours per Battery
Endurance:	Up to 12 km per battery
Overall size:	94 x 65 x 97.5 cm (L x W x H)
Folded size:	65 x 35 x 77.5 cm (L x W x H)

## CH4080 EZEE FOLD DLX

The CH4080 DLX EZee Fold power wheelchair is part of our “top of the line” Deluxe lightweight folding wheelchairs. The CH4080 DLX has 7” front casters with suspension forks, 8” rear wheels and optional elevating leg rests. It is designed for those who want a high quality folding power wheelchair for use primarily indoors and on solid outdoor surfaces.

It’s 18 km range will provide all day independent mobility for even the most active users. Range is dependent on terrain travelled, age and condition of batteries and load weight on wheelchair.



SPECIFICATIONS	
Seat Width:	420mm (16.5")
Seat Height:	510mm (20")
Seat Depth	420mm (16.5")
Backrest Height:	450mm (17.7")
Net Weight:	27.5 kg (60.5 lb)
Weight Capacity:	150kg (330 lb)
Caster:	178mm (7")
Rear Wheel:	203 mm (8")
Battery:	12 v x 2 – 16 AH Lithium
Max speed:	6 km/h (3.75 mph)
Max. Climbing Grade:	12 degrees
Charging period:	8-10 hours per battery
Endurance:	Up to 18 km (11.1847 miles) /battery
Overall size:	1000 x 700 x 890 mm (39.4" x 27.6" x 35")
Folded size:	340 x 700 x 750 mm (13.4" x 27.6" x 29.5")

# G3 POWER WHEELCHAIRS

## CH4085 EZEE FOLD DLX

The CH4085 DLX EZee Fold power wheelchair is part of our “top of the line” Deluxe lightweight folding wheelchairs. The CH4085 DLX has 8” front casters with suspension forks, 12” rear wheels and optional elevating leg rests. It is designed for those who want a high quality folding power wheelchair for use indoors and outdoor surfaces.

It’s 21 km range will provide all day independent mobility for even the most active users. Range is dependent on terrain travelled, age and condition of batteries and load weight on wheelchair.



SPECIFICATIONS	
Seat Width:	458mm (18")
Seat Height:	550mm (21.7")
Seat Depth	420mm (16.5")
Backrest Height:	500mm (19.7")
Net Weight:	33 kg (72.6 lb )
Weight Capacity:	150kg (330 lb)
Caster:	203mm (8")
Rear Wheel:	305mm (12")
Battery:	12 v x 2 – 16 AH Lithium
Max speed:	6 km/h (3 3/4 mph)
Max. Climbing Grade:	12 degrees
Charging period:	8-10 hours per battery
Endurance:	Up to 21 km (13 miles) /battery
Overall size:	1005 x 690 x 990 mm (41.3" x 27.2 x 39")
Folded size:	440 x 690 x 840 mm (17.3" x 27.2 x 33")

## CH4075 (SILVER FRAME) EZEE FOLD

The CH4085 is the largest of our EZee Fold line of portable folding wheelchairs offering a 20" seat and weight capacity of 500 pounds it will take larger adults where they want to go while still providing a folding frame for easy storage and transport.



SPECIFICATIONS	
Seat Width:	508 mm (20")
Seat Height:	540 mm (21.3")
Seat Depth:	508 mm (20")
Backrest Height:	457 mm (18")
Net Weight:	25 kg (55 lb)
Weight Capacity:	227 kg (500 lb)
Caster:	200 mm (8")
Rear Wheel:	305mm (12")
Battery:	24 V, 20 Ah (3.2 kg/7 lb)
Max speed:	8 kph (5 mph)
Max. Climbing Grade:	12 degrees
Charging period:	8-10 hours per battery
Endurance:	Up to 220 km (12.5 miles) /battery
Overall size:	1280 x 656 x 1029 mm (50.4" x 25.8 x 40.5")
Folded Width:	448 x 690 x 840 mm (17.6" x 27.2 x 33")

# G5 POWER WHEELCHAIRS

## CH4070 (SILVER FRAME) EZEE FOLD

Weighing just 37.5 pounds the CH4070 EZee Fold Ultra-Light is the lightest power folding wheelchair on the market we know of. Great for shopping trips and other outings or for use around the house.



SPECIFICATIONS	
Seat Width:	565 mm (22.2")
Seat Height:	900 mm (35.4")
Seat Depth	400 mm (15.75")
Backrest Height:	390 mm (15.4")
Net Weight:	26 kg (57.2 lb.)
Weight Capacity:	100 kg (220 lb)
Caster:	125 mm (5")
Rear Wheel:	200 mm (8")
Battery:	24 volt 10.4 Ah Lithium Ion
Max speed:	4.5 KPH (2.8 MPH)
Static Stability:	≥6°
Charging period:	8-12 hours
Endurance:	Up to 21 km (13 miles) /battery
Overall size:	851 x 559 x 673 mm (33.5" x 22" x 26.5")
Folded size:	406 x 559 x 673 mm (16" x 22" x 26.5")

## EZEE TRAVEL SCOOTER

The Ezee Life EZee Travel Electric scooters are trike style mobility scooters that are super light foldable and easily transportable, and so are the best choice for electric scooters for on-the-move-riders! Folded and stored for the night, it will be ready when you are the next morning. Slide the Ezee Travel into the back of a van with virtually no lifting effort and unfold it again in no time. There are no connectors, latches or mechanisms to open when folding or unfolding. It's that easy!

SPECIFICATIONS	
Drive type:	Rear wheel drive
Maximum Speed:	3.7 mph
Maximum Climbing Slope:	6 degrees (10 %)
Range:	up to 6 miles
Turning Radius:	29"
Ground Clearance:	3.6"
Suspension:	No
Weight capacity:	250 lbs
Battery:	12 A/h 8 kg (18 lbs)
Charging Time:	7 – 10 hours
Rear Wheels:	200 x 50 Solid tires (8" x 2")
Front wheels:	200 x 50 Solid tires (8" x 2")
Overall Length:	39"
Overall Width:	22"
Folding:	Yes
Disassembly:	Partial
Basket:	Underseat tray
Seat Width:	18"
Seat Depth:	16"
Seat Swivel:	No
Lights:	No
Battery gauge:	Yes
Horn:	No
Signals:	No
Rearview Mirror:	No
Overall weight:	72 lbs



# MOBILITY SCOOTER

## EZEE CLASSIC SCOOTER

The Ezee Life EZee Classic Electric scooters are small lightweight mobility scooters that are super light, easily disassemble into four parts or can be folded for easy transportation and provide an excellent choice for electric scooters for on-the-move-riders! Because this mobility scooter easily breaks down into pieces it can be transported in almost any vehicle yet is strong enough to carry adults up to 250 pounds (115 kgs) in weight!

Available in 4 Frame colours: Red, Blue, Black & Silver.

SPECIFICATIONS	
Drive type:	Front wheel drive
Maximum Speed:	4 mph
Maximum Climbing Slope:	6 degrees (10 %)
Range:	up to 6 miles
Turning Radius:	45"
Ground Clearance:	3.6"
Suspension:	No
Weight capacity:	250 lbs
Battery:	12 A/h 8 kg (18 lbs)
Charging Time:	7 – 10 hours
Rear Wheels:	200 x 50 Solid tires (8" x 2")
Front wheels:	200 x 50 Solid tires (8" x 2")
Overall Length:	41"
Overall Width:	22"
Folding:	Yes
Disassembly:	Yes
Basket:	Yes
Seat Width:	18"
Seat Depth:	16"
Seat Swivel:	No
Lights:	No
Battery gauge:	Yes
Horn:	No
Signals:	No
Rearview Mirror:	No
Overall weight:	67 lbs



## EZEE ELITE MOBILITY SCOOTER

The Elite Mobility Scooter is a “Class A” battery powered scooter intended to be used by individuals that are able to walk but suffer from mild mobility limitations. The Elite Mobility Scooter is intended for indoor and outdoor use on pavements or paved footpaths during daylight hours or in well-lit areas.

Available in 4 Frame colours: Red, Blue, Black & Silver.

SPECIFICATIONS	
Drive type:	Front wheel drive
Maximum Speed:	4 mph
Maximum Climbing Slope:	6 degrees (10 %)
Range:	up to 8 miles
Turning Radius:	45"
Ground Clearance:	3.6"
Suspension:	No
Weight capacity:	250 lbs
Battery:	12 A/h 8 kg (18 lbs)
Charging Time:	7 – 10 hours
Rear Wheels:	200 x 50 Solid tires (8" x 2")
Front wheels:	200 x 50 Solid tires (10" x 2")
Overall Length:	41"
Overall Width:	22"
Folding:	Yes
Disassembly:	Yes
Basket:	Yes
Seat Width:	18"
Seat Depth:	16"
Seat Swivel:	No
Lights:	No
Battery gauge:	Yes
Horn:	No
Signals:	No
Rearview Mirror:	No
Overall weight:	59 lbs



# MOBILITY SCOOTER

## EZ5S SCOOTER

The EZ5S mobility scooter is a 3 wheeled portable scooter small enough to be easily loaded into vehicles but powerful enough to be used both indoors and outdoors on solid and relatively smooth surfaces. The EZ5S Scooter disassembles into four pieces without the need of any tools for easy transport and storage with the rear section being the heaviest piece weighing only 31 lbs. (14kgs). With a range of up to 10 miles (16 kms) and a speed of up to 3.5 mph (5.6 kph) the EZ5S mobility scooter will last most users the entire day without the need of recharging. Range and speed will vary depending on user weight, battery condition and terrain.

SPECIFICATIONS	
Drive type:	Rear wheel drive
Maximum Speed:	3.5 mph
Maximum Climbing Slope:	8.00%
Range:	10 miles
Turning Radius:	31.5"
Ground Clearance:	2"
Suspension:	None
Weight capacity:	250 lbs
Battery:	12V 12 AH x 2, Sealed, deep cycle
Charging Time:	10 Hours
Rear Wheels:	7.5" x 2.5" Solid Tires
Front wheels:	7.5" x 2" Solid Tires
Overall Length:	41"
Overall Width:	19.5"
Folding:	no
Disassembly:	Yes
Basket:	Yes
Seat Width:	16.5"
Seat Depth:	14.5"
Seat Swivel:	Yes
Lights:	Yes
Battery gauge:	Yes
Horn:	Yes
Signals:	None
Rearview Mirror:	None
Overall weight:	90.2 lbs



## EZ6S SCOOTER

The EZ6S mobility scooter is a 4 wheeled portable scooter small enough to be easily loaded into vehicles but powerful enough to be used both indoors and outdoors on solid and relatively smooth surfaces. The EZ6S Scooter disassembles into four pieces without the need of any tools for easy transport and storage with the rear section being the heaviest piece weighing only 31 lbs. (14kgs). With a range of up to 10 miles (16 kms) and a speed of up to 3.5 mph (5.6 kph) the EZ6S mobility scooter will last most users the entire day without the need of recharging. Range and speed will vary depending on user weight, battery condition and terrain.

SPECIFICATIONS	
Drive type:	Rear wheel drive
Maximum Speed:	3.5 mph
Maximum Climbing Slope:	8.00%
Range:	10 miles
Turning Radius:	39"
Ground Clearance:	2"
Suspension:	None
Weight capacity:	250 lbs
Battery:	12V 12 AH x 2, Sealed, deep cycle
Charging Time:	10 Hours
Rear Wheels:	7.5" x 2.5" Solid Tires
Front wheels:	7.5" x 2" Solid Tires
Overall Length:	41"
Overall Width:	22"
Folding:	no
Disassembly:	Yes
Basket:	Yes
Seat Width:	16.5"
Seat Depth:	14.5"
Seat Swivel:	Yes
Lights:	Yes
Battery gauge:	Yes
Horn:	Yes
Signals:	None
Rearview Mirror:	None
Overall weight:	99.8 lbs



# MOBILITY SCOOTER

## EZEE ROVER 4 SCOOTER

EZee Life EZee Rover 4 electric scooters are Rugged, Comfortable and Safe. Finally, an electric mobility scooter scooter that is not afraid of the outdoors and can take you places no other scooter can! The Ezee Rover 4 comes apart and easily folds for transport or storage and fits into the trunk of most cars. This Rugged mobility scooter comes with a full light package for optimal safety, and a keyed ignition for added security. A nice wide 18 inch deluxe memory foam seat, and big 10" pneumatic EASY-RIDE tires on the front wheels, and 14" tires on the rear wheels give the Rover 4 scooter a stable comfortable ride on most terrain.

SPECIFICATIONS	
Drive type:	Rear wheel drive
Maximum Speed:	6.8 mph
Maximum Climbing Slope:	8 degrees (10%)
Range:	15.5 miles
Turning Radius:	52"
Ground Clearance:	5"
Suspension:	Yes
Weight capacity:	375 lbs
Battery:	12 Volt 26 Ah Sealed Lead Acid
Charging Time:	7 – 10 hours
Rear Wheels:	14" pneumatic
Front wheels:	10" pneumatic
Overall Length:	54"
Overall Width:	22"
Folding:	Yes
Disassembly:	Yes
Basket:	Yes
Seat Width:	18"
Seat Depth:	18.5"
Seat Swivel:	Yes
Lights:	Yes
Battery gauge:	Yes
Horn:	Yes
Signals:	Yes
Rearview Mirror:	No
Overall weight:	141 lbs




## EZEE AUTO FOLD SCOOTER

The EZee Auto-Fold Scooter is an easy to use power folding portable mobility scooter operated with a remote control or manually. The EZee Auto-Fold has a 300 lb. weight capacity and its standard lightweight Lithium battery provides up to 21.7 km on a single charge (range varies depending on user weight, terrain and battery condition).

The EZee Auto-Fold easily fits into the back of even the smallest hatchback. It is also airline friendly and fully compliant and within regulations for Lithium Ion dry cell batteries.

Available in 4 Frame colours: Red, Orange & Yellow.

SPECIFICATIONS	
Drive type:	Rear wheel drive
Maximum Speed:	6 kph (3.7 mph)
Maximum Climbing Slope:	15 degrees
Range:	13.5 Miles
Turning Radius:	55"
Ground Clearance:	2"
Suspension:	No
Weight capacity:	136 kg (300 lbs.)
Battery:	2 x 24V 6AH lithium-ion
Charging Time:	7 – 10 hours
Rear Wheels:	7" Solid
Front wheels:	6" Solid
Overall Length:	36.9"
Overall Width:	18.1"
Folding:	Yes
Disassembly:	No
Basket:	No
Seat Width:	17"
Seat Depth:	16"
Seat Swivel:	No
Lights:	No
Battery gauge:	Yes
Horn:	Yes
Signals:	No
Rearview Mirror:	No
Overall weight:	53 lbs



## **Benefits of using EZee Life ECVs:**

- 1) Folding and compact, EZee Life Power Wheelchairs & Scooters take up less space than conventional power mobility devices when not in use.
- 2) EZee Life lithium batteries provide all day power and can be charged overnight to be ready for morning use every day.
- 3) Virtually maintenance free with flat-free tires, sealed batteries and designed for long term trouble free use.
- 4) Easy for users to understand and drive.
- 5) Programmable drive controls to restrict speed and power.
- 6) Excellent warranty against manufacturer defects.

## **Benefits of buying your ECVs from Factory Direct Medical:**

- 1) 27 years of medical equipment experience.
- 2) Fast service from our Concord, Ontario warehouse.
- 3) As a medical supply company we are not required to collect sales taxes on mobility devices.
- 4) Registered and licensed as an importer, manufacturer and dealer of medical equipment by Health Canada.
- 5) We have the largest selection of Wheelchair and Scooter ECV models in North America. Get the models most suitable for your facility!





**34 Futurity Gate, Unit # 15 Concord Ontario L4K 1S6**

**wk: 416-739-1267 | Toll Free: 1-888-634-5808 | fx 416-739-7171**

**email: [info@hpurehab.com](mailto:info@hpurehab.com) | [www.hpurehab.com](http://www.hpurehab.com)**