

ATTO ACCESSORIES are created specifically to enhance, support and work with the ATTO mobility scooter in all user modes; while driving, while wheeling along as trolley and when folded for storage or transportation.

**MORE FOR ATTO**

movinglife®


**ATTO**

CARRYALL & CUSHION

Don't think twice, just bring it along!




# HOW TO USE




**1**

Slip the cushion over the seat (use handles to pull if needed)




**2**

Simply pull hook on elastic strap to fold under and connect to ATTO chassis



**3**

One had 'pull to open' Secure funnel access so nothing falls out. Quick 'no zip' flap closing (to ensure nothing falls out zip up when you can)



**4**

Attach to ATTO Padded Cushion dedicated lower zip