



ATTO ACCESSORIES are created specifically to enhance, support and work with the ATTO mobility scooter in all user modes; while driving, while wheeling along as trolley and when folded for storage or transportation.

MORE FOR **ATTO**


movinglife®



ATTO
STANDARD BATTERY


Extra power gets you far!

A | INSTALLING THE BATTERY



1


Insert the battery starting with the back end



2


Push the front down till it clicks into place

B | REMOVING THE BATTERY



1

Press the battery release button



2

Pull the battery out

*** NOTE!** Make sure the battery is properly positioned in place before folding ATTO or switching it on